



RONDELL DISCLOSING PELLETS

Rondell blue disclosing pellets are a convenient and easier way of applying disclosing solution. Just rub the presoaked pellet on the tooth surface and plaque becomes visible.

Plaque that is younger than 3 days colours red and plaque that is older than 3 days is coloured blue.

The solution has a mild peppermint flavour.



***Further information can be obtained from Oraltec NZ Ltd
Freephone 0508 478 222 : Website: www.oraltec.co.nz***



NEW TEMPORARY CEMENT FROM CAVEX

New, improved Cavex temporary cement is a two paste system – a white base paste and yellow catalyst which when mixed together become a light yellow cement.

The cement sets in 4-5 minutes and is a different colour to dentine so it is easy to remove cement traces.



***Further information can be obtained from Oraltec NZ Ltd
Freephone 0508 478 222 : Website: www.oraltec.co.nz***

FENDAWEDGE

Protecting adjacent teeth is easy with new Fenderwedges from Directa Dental. Their unique design combines a plastic interdigital wedge with a built-in vertical stainless steel band that protects adjacent teeth from accidental damage by a bur.

Fenderwedges slip effortlessly into the interproximal space and during use, they ease the teeth apart and aid the insertion of a matrix.

Fenderwedges are available in four sizes and are colour coded for easy identification.



***Further information can be obtained from Oraltec NZ Ltd
Freephone 0508 478 222 : Website: www.oraltec.co.nz***