




## NiTinol Sectional Matrix System

The latest in  
**advanced matrix  
technology** for  
Class II composite  
restorations.



### Consistency

Nitinol ring composition provides improved strength and spring memory retention allowing for more consistent clamping and separation forces across a wide range of tooth sizes.

### Longevity

NiTin's drawn wire Nitinol ring construction, with its superior molecular properties strengthened by a proprietary "Heat Set" process, offers unparalleled longevity versus traditional nickel-titanium.

### Retention

The tip design of the NiTin ring provides significantly improved resistance to "slip off" than other sectional matrix rings.

### Clinical Performance Advantage

Consistently tight contacts with a snug gingival seal creates anatomically correct restorations

////////////////////////////////////

Available exclusively from Oraltec NZ Ltd

**re·invent**  
DENTAL PRODUCTS



## NiTinol Sectional Matrix System

### NiTinol Complete Kit

NiTinol Sectional Matrix System RENTKMN01 :

#### INCLUDES

1	Premolar NiTin Ring	80	Assorted G-Wedge™
1	New Molar NiTin Ring	1	Ring Placement Forceps
80	Assorted Full Curve matrices	1	Band Forceps



### Order Information

	Description	Qty.
<b>Rings</b>		
RENT400	Molar	2
RENT500	Premolar	2

#### Bands


RENTM100	Premolar, 4.4 mm	50
RENTM150	Premolar w/ext, 6.0 mm	25
RENTM175	Molar, 5.6 mm	50
RENTM200	Tall Molar, 6.6 mm	50
RENTM300	Tall Molar w/ext, 8.7 mm	25

#### Wedges

RENTYL	Wedge refill Yellow	100
RENTBL	Wedge refill Blue	100
RENTOR	Wedge refill Orange	100
RENTGR	Wedge refill Green	100
RENTK4	Wedge Kit	400

#### Instruments

RENTMRDF100	Ring Placement Forceps	1
RENTDF68	Band Placement Forceps	1



Ask your Oraltec Representative for a demo or phone 0508-478-222 for further information  
[www.oraltec.co.nz](http://www.oraltec.co.nz) / email: [orders@oraltec.co.nz](mailto:orders@oraltec.co.nz)